



Bioactive glass powder, 0-20 μm

Catlog Number: BGP-10313

• Basic Information

Particle Size	0-20 μm
Chemical Composition	24.5% CaO, 45% SiO ₂ , 6% P ₂ O ₅ , 24.5% Na ₂ O
Color	White
Form	Powder



CD BioSustainable

• Storage

Store in a cool, dry place away from direct sunlight.

• Applications

It can be used as a raw material for the preparation of bioactive glass scaffolds for artificial bone repair materials, as well as for the development of bioactive additives for advanced cosmetics.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA