



Bioactive glass white powder

Catlog Number: BGP-10315

• Basic Information

Particle Size	0-80 μm
Chemical Composition	20% CaO, 53% SiO ₂ , 4% P ₂ O ₅ , 5% MgO, 6% Na ₂ O, 12% K ₂ O
Color	White
Form	Powder



CD BioSustainable

• Storage

Store in a cool, dry place away from direct sunlight.

• Applications

It can be used as a raw material for the preparation of bioactive glass scaffolds for artificial bone repair materials, as well as for the development of bioactive additives for advanced cosmetics.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA