



## Bioactive glass powder, 0.2-20 $\mu\text{m}$

Catlog Number: BGP-10319

### • Basic Information

Particle Size	0.2-20 $\mu\text{m}$
Purity	98.00%
Color	White
Form	Powder



CD BioSustainable

### • Storage

Store in a cool, dry place away from direct sunlight.

### • Applications

It can be used as a high temperature binder for tissue engineering, enamel reconstruction, and bio-coatings.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595  
Email: [info@sustainable-bio.com](mailto:info@sustainable-bio.com)  
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA