



Bioactive glass powder, 50-250 μm

Catlog Number: BGP-10320

• Basic Information

Particle Size	50-250 μm
Purity	98.00%
Color	White
Form	Powder



CD BioSustainable

• Storage

Store in a cool, dry place away from direct sunlight.

• Applications

It can be used as a high temperature binder for tissue engineering, enamel reconstruction, and bio-coatings.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA