



Lecithin

Catlog Number: MDP-11110

• Basic Information

CAS No.	8002-43-5
Purity	≥98.0%
Synonyms	Polyene Phosphatidylcholine
Appearance	Solid
Color	Light brown to yellow
SMILES	<chem>O=P(OCC[N+](C)(C)C)([O-])OC[C@@H](COC([R])=O)OC([R'])=O.[R,R'=fatty acid residues]</chem>
Intended Use	For research use only.



CD BioSustainable

• Package and Storage

Shipping: Room temperature

Storage: Powder: -20°C 3 years; 4°C 2 years

In solvent: -80°C 6 months; -20°C 1 month

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA