



γ -Tocotrienol

Catlog Number: MDP-11335

• Basic Information

CAS No.	14101-61-2
Purity	99.73%
Molecular Weight	410.63
Formula	C ₂₈ H ₄₂ O ₂
Appearance	Liquid (Density: 0.964±0.06 g/cm ³)
Color	Colorless to light yellow
SMILES	<chem>OC1=C(C)C(C)=C2C(CC[C@](C)(C/C=C(C)/CC/C=C(C)/CC/C=C(C)/C)(C)O2)=C1</chem>
Intended Use	For research use only.



CD BioSustainable

• Package and Storage

Shipping: Room temperature

Storage: Pure form -20°C 3 years

In solvent: -80°C 6 months; -20°C 1 month

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA