



Pepstatin Trifluoroacetate

Catlog Number: MDP-11470

• Basic Information

Purity	98.77%
Synonyms	Pepstatin A Trifluoroacetate
Molecular Weight	799.92
Formula	C ₃₆ H ₆₄ F ₃ N ₅ O ₁₁
Appearance	Solid
Color	White to off-white
Intended Use	For research use only.



CD BioSustainable

• Package and Storage

Shipping: Room temperature

Storage: Sealed storage, away from moisture

Powder: -80°C 2 years; -20°C 1 year

In solvent: -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA