



4-Methylpentanoic Acid

Catlog Number: MDP-11734

• Basic Information

CAS No.	646-07-1
Purity	99.56%
Synonyms	Isocaproic Acid
Molecular Weight	116.16
Formula	C ₆ H ₁₂ O ₂
Appearance	Liquid (Density: 0.9392 g/cm ³)
Color	Colorless to light yellow
SMILES	<chem>CC(C)CCC(O)=O</chem>
Intended Use	For research use only.



CD BioSustainable

• Package and Storage

Shipping: Room temperature

Storage: Pure form -20°C 3 years; 4°C 2 years

In solvent: -80°C 6 months; -20°C 1 month

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA