



Ganoderic Acid T-Q

Catlog Number: MDP-11786

• Basic Information

CAS No.	112430-66-7
Purity	99.96%
Synonyms	Ganodermic Acid T-Q
Molecular Weight	510.70
Formula	C32H46O5
Appearance	Solid
Color	White to off-white
SMILES	<chem>CC1(C)C(CC[C@]2(C)C3=CC[C@]4(C)[C@@H]([C@H](C)CC/C=C(C)/C(O)=O)[C@H](OC(C)=O)[C@](C)4C3=CC[C@@]12[H])=O</chem>
Intended Use	For research use only.



CD BioSustainable

• Package and Storage

Shipping: Room temperature

Storage: Powder: -20°C 3 years; 4°C 2 years

In solvent: -80°C 6 months; -20°C 1 month

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA