



Ganoderic Acid D

Catlog Number: PDP-13687

• Basic Information

Appearance	Solid
CAS	108340-60-9
Purity	99.40%
Molecular Weight	514.65
Formula	C ₃₀ H ₄₂ O ₇
Color	White to off-white
SMILES	<chem>CC1(C)C(CC[C@]2(C)C3=C([C@@]4(C(C[C@@H]([C@]4(CC3=O)C)[C@H](C)CC(CC(C)C(O)=O)=O)=O)C)[C@@H](O)C[C@@]12[H])=O</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

• Package and Storage

Storage: Powder: -20°C 3 years; 4°C 2 years

In solvent: -80°C 6 months; -20°C 1 month

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA