



Soyasapogenol A

Catlog Number: PDP-15100

• Basic Information

Appearance	Solid
CAS	508-01-0
Purity	99.06%
Molecular Weight	474.72
Formula	C30H50O4
Color	White to off-white
SMILES	<chem>C[C@]12[C@]3(C([C@@]4([H])[C@](C)([C@@H]([C@H](O)C(C)(C)C4)O)CC3)=CC[C@]1([H])[C@@]5([C@@]([C@](C)([C@@H](O)CC5)CO)([H])CC2)C)C</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

• Package and Storage

Storage: 4°C, protect from light

In solvent: -80°C, 6 months; -20°C, 1 month (protect from light)

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA