



Cyanidin 3-sambubioside Chloride

Catlog Number: PDP-15130

• Basic Information

Synonyms	Cyanidin-3-O-sambubioside Chloride
Appearance	Solid
CAS	33012-73-6
Purity	98.50%
Molecular Weight	616.95
Formula	C ₂₆ H ₂₉ ClO ₁₅
Color	Brown to black
SMILES	<chem>OC1=CC(O)=CC2=C1C=C(C(C3=CC(O)=C(C=C3)O)=[O+])2O[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)[C@H]5[C@@H]([C@H]([C@@H](CO5)O)O)O.[Cl-]</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

• Package and Storage

Storage: 4°C, sealed storage, away from moisture and light

In solvent: -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture and light)

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA