



Soyasapogenol C

Catlog Number: PDP-15361

• Basic Information

Appearance	Solid
CAS	595-14-2
Purity	≥99.0%
Molecular Weight	440.70
Formula	C30H48O2
Color	White to off-white
SMILES	<chem>C[C@]12[C@]3(C([C@@]4([H])[C@@](CC3)(C=CC(C)(C4)C)C)=CC[C@]1([H])[C@@]5([C@@]([C@](C)([C@H](CC5)O)CO)([H])CC2)C)C</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

• Package and Storage

Storage: 4°C, protect from light

In solvent: -80°C, 6 months; -20°C, 1 month (protect from light)

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA