



Tokinolide B

Catlog Number: PDP-15537

• Basic Information

Appearance	Solid
CAS	112966-16-2
Purity	98.44%
Molecular Weight	380.48
Formula	C ₂₄ H ₂₈ O ₄
Color	White to off-white
SMILES	<chem>O=C(O1)C2=C([C@]31[C@@]4(CC5)C(C(O/C4=C\CC C)=O)=C[C@]5([H])[C@H]3CCC)CCC=C2</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

• Package and Storage

Storage: -20°C, protect from light

In solvent: -80°C, 6 months; -20°C, 1 month (protect from light)

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA