



Ginsenoside Ra3

Catlog Number: PDP-15605

• Basic Information

Appearance	Solid
CAS	90985-77-6
Purity	98.24%
Molecular Weight	1241.41
Formula	C ₅₉ H ₁₀₀ O ₂₇
Color	White to off-white
SMILES	<chem>CC12C3(C(C(O)CC1C4(C(C(C)(C(OC5C(C(C(O)C(CO)O5)O)OC6OC(C(O)C(O)C6O)CO)CC4)C)CC2)C)C(C(C/C=C(C)/C)(C)OC7OC(C(O)C(O)C7O)COC8C(C(C(O)C(CO)O8)OC9C(C(C(O)CO9)O)O)O)CC3)C</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

• Package and Storage

Storage: -20°C, protect from light

In solvent: -80°C, 6 months; -20°C, 1 month (protect from light)

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA