



Tanshinone IIB

Catlog Number: PDP-15881

• Basic Information

Appearance	Solid
CAS	17397-93-2
Purity	≥99.0%
Molecular Weight	310.34
Formula	C ₁₉ H ₁₈ O ₄
Color	White to off-white
SMILES	<chem>O=C(C1=C2C=CC3=C1CCC[C@@]3(CO)C)C(C4=C2OC=C4C)=O</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

• Package and Storage

Storage: 4°C, protect from light

In solvent: -80°C, 6 months; -20°C, 1 month (protect from light)

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA