



Saikogenin A

Catlog Number: PDP-16154

• Basic Information

Appearance	Solid
CAS	5092-09-1
Purity	99.99%
Molecular Weight	472.70
Formula	C30H48O4
Color	White to off-white
SMILES	<chem>C[C@]12[C@]3(C(C=C[C@]1([H])[C@@]4([C@@]([C@@](C)([C@@H](O)CC4)CO)([H])CC2)C)=C5[C@](C O)(CCC(C)(C)C5)[C@@H](O)C3)C</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

• Package and Storage

Storage: Powder: -20°C 3 years; 4°C 2 years

In solvent: -80°C 6 months; -20°C 1 month

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA