



## Sinigrin Hydrate

Catlog Number: PDP-16387

### • Basic Information

Appearance	Solid
CAS	64550-88-5
Purity	99.77%
Molecular Weight	415.48
Formula	C <sub>10</sub> H <sub>18</sub> KNO <sub>10</sub> S <sub>2</sub>
Color	White to off-white
SMILES	<chem>O[C@@H]1[C@@H](O)[C@H](S/C(CC=C)=N/OS(=O)(O[K])=O)[C@H](CO)[C@H]1O.O</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

### • Package and Storage

Storage: -20°C, sealed storage, away from moisture

In solvent: -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595  
Email: [info@sustainable-bio.com](mailto:info@sustainable-bio.com)  
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA