



## Ganoderic Acid C1

Catlog Number: PDP-16502

### • Basic Information

Appearance	Solid
CAS	95311-97-0
Molecular Weight	514.65
Formula	C <sub>30</sub> H <sub>42</sub> O <sub>7</sub>
Color	White to light yellow
SMILES	<chem>C[C@]12C3=C(C(C[C@@]1([C@@](CC2=O)([H])[C@H](C)CC(C[C@@H](C)C(O)=O)=O)C)=O)[C@@]4([C@@](C(C)(C(CC4)=O)C)([H])C[C@@H]3O)C</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

### • Package and Storage

Storage: 4°C, protect from light

In solvent: -80°C, 6 months; -20°C, 1 month (protect from light)

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595  
Email: [info@sustainable-bio.com](mailto:info@sustainable-bio.com)  
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA