



3,3'-Bisdemethylpinoresinol

Catlog Number: PDP-16565

• Basic Information

Appearance	Solid
CAS	340167-81-9
Purity	98.96%
Molecular Weight	330.33
Formula	C18H18O6
Color	White to off-white
SMILES	<chem>OC1=C(O)C=CC([C@@H]2[C@]3([H])[C@@](CO2)([H])[C@@H](C4=CC(O)=C(C=C4)O)OC3)=C1</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

• Package and Storage

Storage: 4°C, sealed storage, away from moisture

In solvent: -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA