



Pinolenic Acid

Catlog Number: PDP-16655

• Basic Information

Synonyms	PNLA
Appearance	Liquid
CAS	16833-54-8
Purity	≥98.0%
Molecular Weight	278.43
Formula	C ₁₈ H ₃₀ O ₂
Color	Colorless to light yellow
SMILES	<chem>CCCCC/C=C\C/C=C\C/C=C\C/C=C\CCCC(O)=O</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

• Package and Storage

Storage: Solution, -20°C, 2 years

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA