



Tetramethylkaempferol

Catlog Number: PDP-17084

• Basic Information

Appearance	Solid
CAS	16692-52-7
Purity	99.93%
Molecular Weight	342.34
Formula	C ₁₉ H ₁₈ O ₆
Color	White to off-white
SMILES	<chem>O=C1C(OC)=C(C2=CC=C(OC)C=C2)OC3=CC(OC)=CC(OC)=C13</chem>
Intended Use	For research and further manufacturing use only.



CD BioSustainable

• Package and Storage

Storage: -20°C, protect from light

In solvent: -80°C, 6 months; -20°C, 1 month (protect from light)

Shipping: Room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA