



Polyethylenimine Hydrochloride, Linear (MW 4,000)

Catlog Number: BBP-0598

• Basic Info

CAS Number: 49553-93-7

Molecular Weight: 4,000 (~2,200 free base)

Color: White to off-white

Physical State: Solid



CD BioSustainable

• Properties

Appearance: White to off-white free flowing solid

Solubility: Cold and room temperature water Insoluble in: Common organic solvents (ethanol, acetone, tetrahydrofuran)

• Transport and Storage

Handling: Glove, chemical goggles & fume hood

Storage: Store at room temperature.

• Security Information

Hazards: Irritant

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA