



Pilinia Earleae

Catlog Number: AB-6403

• Basic Info

Class: Chlorophyceae

Strain: Pilinia earleae

Medium: Enriched Seawater Medium

Type Culture: Yes

Synonyms: Pilinia sp.



CD BioSustainable

• Culture Condition

Temperature: 20°C

Light Source: Mix of warm-white & cool-white fluorescent lamps.

Intensity: 3200 lux

Periodicity: 12/12 h L/D

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA