



Poly(acrylic acid), 25% soln. in water [Mw ~345,000]

Catlog Number: BBP-0679

• Basic Info

CAS Number: 9003-01-4

Molecular Weight: ~345,000 g/mol



CD BioSustainable

• Properties

Refractive Index: n_{25/D} 1.527

Viscosity: 400-1200 cps

• Transport and Storage

Handling: Exercise normal care

Storage: Store at room temperature.

• Security Information

Hazards: Harmless-use normal precautions

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA