



Polyacrylamide, Mw 10,000, 50 wt. % in H₂O

Catlog Number: BBP-0695

• Basic Info

CAS Number: 9003-05-8

Molecular Weight: 10,000 g/mol

Physical State: Liquid



CD BioSustainable

• Properties

Appearance: Viscous liquid

Viscosity: 15% AQ 5-20 cps

Solubility: Water, morpholine

Glass Transition Temperature: 165°C

• Transport and Storage

Handling: Gloves & chemical goggles

Storage: Store at room temperature.

• Security Information

Hazards: Irritant

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA