



## Wheat genotyping chip

Catlog Number: PAAGM-8148

### • Basic Information

Species	Wheat
Specification	96 Array Plate
Bitmap	817 K
Expiration Date	1 Year



CD BioSustainable

### • Application

It can be used for complex trait studies such as genetic diversity analysis, association mapping, genome-wide analysis and selection.

### • Storage

Store at -20°C~-80°C.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595  
Email: [info@sustainable-bio.com](mailto:info@sustainable-bio.com)  
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA