



Total Nitrogen Special Reagent (Economical Type)

Catlog Number: WQMR-8232

• Overview

Features: The reagents are easy to prepare. Dangerous chemical reagents will be used and require professionals to operate.



CD BioSustainable

• Basic Information

| | |
|-----------------|-----------|
| Size | 100 times |
| Detection Range | 0~25 mg/L |

• Package and Storage

Shelf Life: 6 months

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA