



COD Detection Prefabricated Reagent (Low Range)

Catlog Number: WQMR-8245

• Overview

Features: Safety: Each sample is paired with a single reagent, eliminating the need for complex preparation and reducing contact with hazardous chemicals. The entire process is automated, ensuring personnel safety.

Simplicity: The testing process is straightforward and easy to operate, requiring no specialized training, making it simple for anyone to perform.

Efficiency: With an easy-to-use, ready-to-go operation mode, sample testing can be completed in just a few minutes.

Environmental friendliness: The reagent usage is minimal, generating only a fraction of the waste liquid compared to traditional methods.

Reliability: Strict quality control and advanced production processes ensure reagent consistency and testing accuracy.

High efficiency: Testing efficiency is greatly improved, enabling each person to complete hundreds of samples per day.

Cost control: Reduces complex procedures, boosts lab efficiency, minimizes waste liquid discharge, saves time, reduces labor, lowers risks, and decreases emissions, ultimately lowering overall costs.

Strong compatibility: The reagent is packaged in high borosilicate glass, making it compatible with various domestic and imported multi-parameter water quality testing devices.

Application: Applicable to the determination of chemical oxygen demand (COD) in surface water, groundwater, domestic sewage and industrial wastewater



CD BioSustainable

• Basic Information

Size	25/150 items
Detection Range	10 ~ 150 mg/L
Determination Time	20 min

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA



COD Detection Prefabricated Reagent (Low Range)

Catlog Number: WQMR-8245

Anti-chloride Ion 1000 mg/L
Interference

• Package and Storage

Shelf Life: 12 months at room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA