



## Total Phosphorus Detection Prefabricated Reagent (High Range)

Catlog Number: WQMR-8250

### • Overview

Application: Applicable to the determination of total phosphorus in surface water, groundwater, domestic sewage and industrial wastewater



CD BioSustainable

### • Basic Information

Size	25/150 items
Detection Range	0.2 ~ 30 mg/L
Type	Pre-assembled reagent.
Determination Time	45 min

### • Package and Storage

Shelf Life: Store at low temperature for 12 months.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595  
Email: [info@sustainable-bio.com](mailto:info@sustainable-bio.com)  
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA