



COD Detection Prefabricated Reagent (Ultra-low Range)

Catlog Number: WQMR-8253

• Overview

Application: Applicable to the determination of chemical oxygen demand (COD) in surface water, groundwater, domestic sewage and industrial wastewater



CD BioSustainable

• Basic Information

Size	25/150 items
Detection Range	0 ~ 40 mg/L
Type	Pre-assembled reagent.
Determination Time	20 min
Anti-chloride Ion	1000 mg/L
Interference	

• Package and Storage

Shelf Life: 12 months at room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA