



## Ferrous Acid Detection Prefabricated Reagents

Catlog Number: WQMR-8266

### • Basic Information

Size	25 times
Detection Range	0~10 mg/L
Type	Pre-assembled reagent.



CD BioSustainable

### • Package and Storage

Shelf Life: 12 months at room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595  
Email: [info@sustainable-bio.com](mailto:info@sustainable-bio.com)  
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA