



Silicon Detection Reagents (High Range)

Catlog Number: WQMR-8269

• Basic Information

Size	100 times
Detection Range	1~80 mg/L
Type	Pre-assembled reagent.



CD BioSustainable

• Package and Storage

Shelf Life: 12 months at room temperature.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA