



## Total Nitrogen Test Reagents

Catlog Number: WQMR-8289

### • Overview

Description: The operation is simple, and the determination principle aligns with methods recognized by the environmental protection industry.

Features: Safe and reliable

Standardized packaging, in line with requirements.

Easy to operate

Semi-finished reagents, easy to operate.

Good quality and low price

The reagents are economical and the testing cost is greatly reduced.

High stability

With hundreds of thousands of test data as reference, the measurement range is wide.

Application: It is widely used in sewage discharge of various industries, domestic sewage discharge, laboratories, etc.



CD BioSustainable

### • Basic Information

Size 100 times

Detection Range 0~100 mg/L

### • Package and Storage

Storage: Store at 2-8°C.

Shelf Life: 6 months.

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595  
Email: [info@sustainable-bio.com](mailto:info@sustainable-bio.com)  
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA