



Surrogate Mix 2, 1000 µg/mL

Catlog Number: WQMR-9258

• Specifications

Size: 1 mL

Property: Colourless Liquid

Purity: 2000 µg/mL in Methanol

Product Tags: Standard solution

Intended Use: For research use only!



CD BioSustainable

• Package and Storage

Storage: -20°C±5°C

Shelf Life: More than 1 year

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA