



Cyanide Solution

Catlog Number: WQMR-9354

• Specifications

Size: 50 mL

Purity: 50 mg/L in 0.1 mol/L NaOH

Product Tags: Standard solution

Intended Use: For research use only!



CD BioSustainable

• Package and Storage

Storage: 2°C-8°C

Shelf Life: More than 1 year

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595
Email: info@sustainable-bio.com
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA