



## Gonyautoxin 1 and Gonyautoxin 4

Catlog Number: WQMR-9367

### • Specifications

---

Size: 0.5 mL

Property: Colourless Liquid

Purity: 5-25 µg/mL in Hydrochloric Acid & Water

Product Tags: Standard solution

Intended Use: For research use only!



CD BioSustainable

### • Package and Storage

---

Storage: 2°C-8°C

Shelf Life: More than 1 year

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595  
Email: [info@sustainable-bio.com](mailto:info@sustainable-bio.com)  
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA