



## Hafnium oxide $\text{HfO}_2$

Catlog Number: NO-9953

### • Specifications

Size: 1  $\mu\text{m}$  100 g; 1  $\mu\text{m}$  500 g; 1  $\mu\text{m}$  1000 g

Form: Powder

Melting Point: 2700°C

Dielectric Constant: 25

Average Particle Size: 1  $\mu\text{m}$

Purity: 99%

Specific Surface Area: 50  $\text{m}^2/\text{g}$

Bulk Density: 8  $\text{g}/\text{cm}^3$

Crystaltype: Quadripartite

Color: White



CD BioSustainable

### • Applications

This product is used as a catalyst support material in catalytic reactions.

### • Product Features

Excellent electrical insulation properties

High melting point

Chemical stability.

High dielectric constant

Good thermal stability

CD BioSustainable is a reputable and professional company in the industry. We provide high-level technical services and total solutions to deliver only the most satisfactory results for our clients.

Tel: 1-631-637-0420; 1-631-533-0595  
Email: [info@sustainable-bio.com](mailto:info@sustainable-bio.com)  
SUITE 206, 17 Ramsey Road, Shirley, NY 11967, USA